

The ACS NCP-2™ External Counterpulsation (ECP) Device



Therapeutic Augmentation is as much about the administration of therapy as it is technology. **ACS** is the first to introduce a comprehensive External Counterpulsation (ECP) Treatment Report. This treatment report enables physicians to quickly **verify** and assess patient progress throughout the 35-day treatment session. Various trends and graphs allow physicians to track and monitor staff interaction. Most importantly, physicians now have a method for **verifying** that ECP treatment is properly administered.

The **ACS NCP-2** maintains a user-friendly interface and incorporates new, innovative reporting capabilities. In addition to a comprehensive ECP Treatment Report, the **ACS NCP-2**

enables staff to monitor patients while utilizing various treatment screens. Staff can now manage time more effectively by performing other tasks. The **ACS NCP-2** is the solution for optimizing patient care and increasing productivity!

Benefits of the **ACS NCP-2™** include:

- ♥ **Automatic Treatment Report** – ACS is the first to generate an automated ECP Treatment Report. Physicians can quickly verify the administration of therapy and assess progress to optimize patient care.
- ♥ **Windows XP™ Operating System** – The ACS NCP-2 operates on Windows XP™ platform enabling a variety of networking, Internet, and communications capabilities.
- ♥ **Dual Utility Platform** – Observe ECP treatment, download and analyze Holter Monitoring reports, or review ICG hemodynamic data on the same operating system. Multi-tasking made easy!
- ♥ **ACS Remote Assistance Program** – The NCP-2 can provide immediate on-line technical or educational support with ACS while in operation.
- ♥ **Mobile Control Console** – This control station is designed with mobility to accommodate most treatment rooms.
- ♥ **Large Format LCD Display** – The ACS Treatment Report is displayed on a 17" LCD screen for ease of operation.
- ♥ **HP Laser Printer** – Utilize the HP Laser printer to print high resolution Full Disclosure and the ECP Treatment Report.
- ♥ **Adjustable Treatment Platform** – Operators can adjust platform levels for patient comfort and compliance. The raised edges of the treatment platform provide increased stability for CHF and angina patients.
- ♥ **Retractable Hoses** – Contour to the patient's body for patients of varying heights.
- ♥ **Interchangeable Cuffs** – Made of a durable ballistic nylon. The improved design includes additional quilting to increase comfort. Optimize therapy by custom fitting cuffs to each patient. There are three sizes available which can be used interchangeably.